## GENTLENEST USER GUIDE

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## EN

Dear customer!

Thank you for choosing a product by Kinderkraft.

Our designs are child-oriented – safety and quality always comes first, stay comfortable knowing that you have made the best choice.

**PLEASE READ THE USER GUIDE CAREFULLY BEFORE USING THE DEVICE.**

### WARNINGS

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

Use only an external supply with the following specifications: 5,0 V  $\overline{\text{---}}$ ; 1,0 A, with USB type C plug.

Do not expose a battery pack or appliance to fire or excessive temperature.

Metal objects can cause short circuits if they come into contact with the battery terminals.

This appliance contains batteries that are non-replaceable. When the battery is at end of life, the appliance shall be properly disposed of.

Do not put the breast pump main unit into water or other solutions to clean (or into the sterilizer disinfection).

Can not be used when charging.

The device is intended for household use only and is not suitable for commercial use.

Technical data:

Operating temperature range: 0°C - 40°C, Charging temperature range: 20°C - 30°C

### 1. List of Parts (Fig. A):

Nr	Part Name	Quantity
1	Lid	2
2	Membrane	8
3	Milk Collector	2
4	Cover	2
5	Check Valve	8
6	Hose Connection Port	2

7	Main Unit	1
8	Hose Connector	1
9	Silicone Hoses	3
10	Silicone Flange 21 mm	2
11	Silicone Flange 19 mm	2
12	Silicone Flange 17 mm	2
13	Silicone Flange 24 mm	2
14	Storage Bag for Breast Pump	1

## 2. Assembly:

Install the membrane (2) onto the cover (4) (Fig. B) and press it firmly into place until fully seated.

Install the check valve (5) into the cover (4) and ensure that it is correctly positioned (Fig. C).

Next, connect the milk collector (3) to the cover (4) and press it firmly with both hands to ensure that the components do not separate (Fig. E). Insert the long silicone tube (9) and the short silicone tubes (9) into the appropriate positions of the connection port (6). Connect the long silicone tube (9) to the tube connection port (6) and the tube connector (8), and connect the short silicone tubes (9) to the tube connector (8) and the lids (1) (Fig. F).

**Double use:** Hold the tube connection port (6) and insert it into the main unit (7). Then insert two lids (1) into the appropriate positions of the milk collectors (3) (Fig. H). Ensure that all components are correctly installed.

**Single use:** Hold the tube connection port (6) and insert it into the main unit (7). Then insert one lid (1) into the appropriate position of the milk collector (3), and insert the other lid into the tube connector (8) (Fig. G). Ensure that all components are correctly installed.

## 3. Charging the Breast Pump:

Connect the Type-C charging cable to the charging port on the top of the main unit (7), then plug the power adapter (5 V / 1 A) into a power outlet (Fig. J).

The display remains on during charging. After charging is complete, the device will enter sleep mode. Once the power adapter is disconnected, the device will turn off.

- One full charge allows approximately 5–6 cycles. One operating cycle lasts about 30 minutes.
- When the battery level is low, the display will flash until the device automatically shuts down when the battery is depleted.
- Charging time is approximately 2 hours.

**Warning:** The device must not be used while charging.

## 4. Control Panel Description (Fig. K):

1	Suction Level Indicator	6	Increase Suction Level
2	Pumping Mode Indicator	7	Battery Level Indicator
3	Mixed Mode Indicator	8	Decrease Suction Level
4	Stimulation Mode Indicator	9	Mode Selection / Change Button
5	Power On/Off Button		

## Main Unit Operation:

- Press and hold the power button for 2 seconds to turn on the breast pump.
- At first use, the device is set by default to stimulation mode at level 3. The display is on.
- A short press of the power button pauses the current operation; the display will start flashing.
- Another short press resumes operation; the display remains lit.
- After 20 seconds of inactivity, the display will turn off and the device will enter sleep mode.
- Press any button to wake the device.
- Press and hold the power button for 2 seconds to turn off the device.

- If no operation is performed within 5 minutes while paused, the breast pump will automatically shut down.

### 5. Breast Pump Modes (Press the mode selection button to switch between functions):

	Stimulation Mode: suction strength adjustable from level 1 to level 9
	Mixed Mode: suction strength adjustable from level 1 to level 9
	Pumping Mode: suction strength adjustable from level 1 to level 9

### 6. Usage Instructions:

- The size of the flange does not affect the amount of milk expressed, but it does affect user comfort.
- This product is supplied with silicone flanges in four sizes (17, 19, 21, and 24 mm). Make sure to use the most comfortable fit.
- The flange should be centered over the nipple (Fig. I). Ensure the selected flange fits properly.
- Place the device inside a nursing bra, put on the bra, hang the main unit around your neck, and begin pumping milk (Fig. I). The device can be used while sitting or standing.
- If pain or discomfort occurs, turn off the device and adjust the position of the flange (Fig. J). Check that there is no gap between the breast and the silicone flange. Restart the device at the lowest suction level. If pain or discomfort persists, stop using the product immediately and consult a physician.

### 7. Suction Level Adjustment:

During operation, the suction level can be adjusted at any time by pressing the “+” or “-” buttons. Each press increases or decreases the level by one.

- After pressing “+” or “-”, the main unit will temporarily pause operation.
- Adjusting the suction level during operation does not affect the total 30-minute operating cycle.
- Press and hold the “+” or “-” button to change the suction level quickly.
- The product features a memory function: the default mode will be set to the last mode used.

### 8. Pouring the Milk:

After pumping, remove the breast pump and point the flange upward to prevent milk from spilling. Prepare a storage bottle and pour the milk according to the correct pouring method shown in Fig. M.

### 9. Storage of Breast Milk:

Temperature – range	Storage time for a healthy baby	Storage time for a sick or premature baby
Temperature up to 15°C	Up to 24 hours	Up to 1 hour
Cooling bag with ice packs – for transport only	Up to 24 hours	Up to 24 hours
	2–5 days	96 hours
Refrigerator (back wall) (+4°C) – do not store in the refrigerator door	2 weeks	1 week
Refrigerator freezer compartment (shared door of refrigerator and	3–4 months	Several weeks

freezer) (-10°C)		
Fridge-freezer (separate doors for refrigerator and freezer) (-14°C)	6 months or longer	3 months

- Because breast milk expands when frozen, it is recommended to fill bottles only 3/4 full before freezing.
- It is recommended that each stored portion be approximately 180 ml – this allows faster thawing and prevents waste.
- Label each storage container with the date.
- Thaw breast milk for several hours in the refrigerator or in warm water at about 40°C. When thawing in warm water, gently swirl the milk to prevent fat separation.
- Thawed breast milk must not be refrozen and must be consumed within 24 hours.


Note:

- Do not store breast milk in the refrigerator door.
- To preserve nutrients, frozen milk should be thawed in the refrigerator the night before use.
- Do not heat frozen milk in a microwave, as this may cause nutrient loss and risk burning the baby.
- Do not mix fresh breast milk with thawed breast milk.

### 10. Cleaning and Disinfection (fig. N):

Keep the product clean. Do not wash or disinfect the charging cable or the main unit in water; wipe the surface with a paper towel. If the device is accidentally submerged in water or comes into contact with another liquid, stop using it immediately and contact customer service. Wash all other breast pump components (fig. O) only by hand using a mild detergent: Before first use and after every use, rinse under running water, wash with mild soap, and rinse thoroughly with warm water. These components may be disinfected with steam or in boiling water; after the water reaches boiling, turn off the heat and wait three minutes before disinfection. Allow all parts to dry completely before reassembly. It is recommended to replace the membrane and the check valve every 3 months.

### DISPOSAL OF WASTE EQUIPMENT AND BATTERIES

The symbol of the crossed-out wheeled bin placed on equipment, batteries, packaging or accompanying documents indicates that the product must not be disposed of with other waste. It is the user's responsibility to hand over waste equipment and batteries to a designated collection point for proper treatment. Information about the available waste electrical equipment and battery collection system can be found in the shop's information point and in the municipal/commune office.

Each household is a user of electrical and electronic equipment and batteries and therefore a potential producer of waste hazardous to

humans and the environment due to the presence of hazardous substances, mixtures and components in the equipment and batteries. On the other hand, waste equipment is a valuable material from which we can recover raw materials such as copper, tin, glass, iron and others. Proper handling of waste equipment and batteries prevents negative consequences for the environment and human health!

### 11. WARRANTY

The full text of the Terms and Conditions of the Warranty is available on the website [WWW.KINDERKRAFT.COM](http://WWW.KINDERKRAFT.COM)

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